
## 5<sup>th</sup> International Meeting on Materials Science for Energy Related Applications

September 25-26, 2025, Belgrade



## Automatizing data storage, analysis, and sharing for the electrochemical $CO_2$ reduction – the FAIR case of SuPERCO2


Dr. Federico Dattila

Department of Applied Science and Technology

Polytechnic of Turin



General meeting – Prof. Simelys Hernández's group 17/10/25, Turin



# The FAIR problem of CO<sub>2</sub>R


## Publications on CO<sub>2</sub>R

#### POLAROGRAPHIC REDUCTION OF CARBON DIOXIDE

Sir:

We have been engaged for some time in the study of the reduction of carbon dioxide at the dropping mercury cathode. On the basis of a large number of polarograms obtained by means of a Heyrovsky Polarograph of the Sargent Co. with solutions of carbon dioxide in 0.1 molar tetramethylammonium chloride, we have arrived at the following definite conclusions: carbon dioxide exhibits well-defined reduction waves with a half-wave potential of remarkable constancy:  $-2.24 \pm 0.01$  volt referred to the saturated calomel electrode. This is the average of ten independent determinations with varying amounts of carbon dioxide. When hydrogen is passed through the substituted ammonium salt solution before the dissolution of carbon dioxide, the waves are equally well defined, but the half-wave potential is somewhat less negative:  $-2.18 \pm 0.02$ volt referred to the saturated calomel electrode. This is the average of seven independent determinations, with varying amounts of hydrogen and carbon dioxide.


Van Rysselberghe, P.; Alkire, G. J. J. Am. Chem. Soc. **1944**, 66, 1801.






## The need for FAIR protocols

Findability / Accessibility / Interoperability / Reusability









### The limit of current repositories

### **Keywords:**

### **Keywords:**


Electrochemical CO<sub>2</sub> reduction Electrochemical CO<sub>2</sub> reduction, Zenodo










### Al-based solution: MERMaid





## The SuPERCO2 concept







## Synthesis

Repository: <a href="mailto:github.com/FedeDat/SuPERCO2">github.com/FedeDat/SuPERCO2</a>

item: data/CuO/synthesis

### Internal logbook

**Precursors** 

**Protocol** 

1	Precursor	CAS number	C (M)	C (mol)	V (mL)	m (g)	Molar mass (g/mol)
2	CuCl2 2H2O	10125-13-0	0.02	0.00	250.00	0.60	134.45
3	NaOH	1310-73-2	6.08	0.30	50.00	12.17	40.00
4	milliQ		55.50	13.88	250.00	249.96	18.02

#### **Equipment**


- · A Peristaltic Pump was used to control the flow rate during addition of the sodium hydroxide solutions.
- · Magnetic stirrers were essential for ensuring uniform mixing of the reactants, thus they were employed whenever a solution was involved.
- A centrifuge was necessary for separating the precipitated particles from the reaction mixture

#### Synthesis procedure

Synthesis carried out on 03/07/2024.

The synthesis of CuO was a simple water-based precipitation, followed by calcination (Table 1). 250 mL of milliQ water were added to 600 mg of copper chloride (CuCl2···2H2O) under vigorous magnetic stirring to achieve a 20 mM solution. Then a solution of 6 M NaOH was formed by mixing 12.1674 g of NaOH into a 50 mL volumetric flask. After this, half the NaOH solution (25 mL) was added to the copper chloride solution via the peristaltic pump with a fixed flow rate of 25 mL per minute. As soon as NaOH was inserted, a light blue copper hydroxide precipitate appeared (Figure 7). The precipitate was centrifuged at 5000 rpm and cleaned with milliQ water. Four additional cleaning/centrifugation cycles were carried out, where the solvent was promptly removed. Then the resulting precipitate was transferred to a standard oven at 200 °C for 2 hours where the hydroxide decomposed to CuO.






## Characterization


Repository: <a href="mailto:github.com/FedeDat/SuPERCO2">github.com/FedeDat/SuPERCO2</a>

item: data/24-10-30\_CuO/characterization

## Scanning Electron Microscopy

**Post-processing** 








## **Energy Dispersive X-ray**

### **Post-processing**



```
def cm_to_inch(cm):
    inch=float(cm/2.54)
    return inch
def AgCl_to_RHE(V_AgCl,pH,ref):
    if ref == "Ag/AgCl":
        V RHE=+0.197+V AgCl+0.059*pH
        return V RHE
# Read the content of the file
url_EDX = "https://raw.githubusercontent.com/FedeDat/SuPERCO2/main/data/24-10-30_CuO/characterization/EDX.csv"
data_EDX = pd.read_csv(url_EDX, usecols=lambda column: pd.notnull(column))
data_EDX.dropna(axis=1, how='all', inplace=True)
#Calculate ratio between metal and oxygen
data EDX[''+str(data EDX.columns[0])+'/'+str(data EDX.columns[1])+'']=data EDX[data EDX.columns[0]]/data EDX[data EDX.columns[0]]
data_EDX.index=list(range(1, len(data_EDX)+1))
ax1 = data_EDX.iloc[:,:-1].plot(kind='bar', stacked=True, figsize=(cm_to_inch(10), cm_to_inch(8)))
plt.legend(loc='upper center')
ax2 = ax1.twinx() # Create a second y-axis
ax2.plot(data EDX.index-1, data_EDX[data EDX.columns[-1]], linestyle='None', marker='o', markersize=8, markerfacecolor='red'
# Adding titles and labels
ax1.set_xticklabels(list(range(1, len(data_EDX)+1)), rotation=0, ha='right')
```


### X-ray Diffraction

### **Post-processing**

#### SuPERCO2 / data / 24-10-30 CuO / characterization / XRD.csv **Preview** Code Blame 2694 lines (2694 loc) · 46 KB Angle Intensity 20.01513028 862 848 20.04139085 20.06765141 881 856 20.09391198 853 20.12017254 891 20.14643311 34 20.17269367 876

```
def AgCl_to_RHE(V_AgCl,pH,ref):
   if ref == "Ag/AgC1":
        V_RHE=+0.197+V_AgCl+0.059*pH
        return V RHE
url XRD = "https://raw.githubusercontent.com/FedeDat/SuPERCO2/main/data/24-10-30 CuO/characterization/XRD.csv"
data_XRD = pd.read_csv(url_XRD, header=None, usecols=[0, 1], names=['Angle', 'Intensity'])
k_XRD = float(data_XRD[data_XRD.apply(lambda row: row.astype(str).str.contains('K-Alpha1 wavelength')).any(axis=1)]['Intensit
n XRD = int(data XRD[data XRD.apply(lambda row: row.astype(str).str.contains('No. of points')).any(axis=1)]['Intensity'])
data XRD sl = data XRD.iloc[(data XRD[data XRD.apply(lambda row: row.astype(str).str.contains('Angle')).any(axis=1)].index[0]
data_XRD_sl = data_XRD_sl.astype(float)
data_XRD_sl['Intensity'] = data_XRD_sl['Intensity'] - np.min(data_XRD_sl['Intensity'])
data XRD sl['Intensity']=data XRD sl['Intensity']/np.max(data XRD sl['Intensity'])
data XRD sl['Intensity savgol'] = savgol filter(data XRD sl['Intensity'], window length=10, polyorder=1)
k_XRD = float(data_XRD[data_XRD.apply(lambda row: row.astype(str).str.contains('K-Alpha1 wavelength')).any(axis=1)]['Intensit
n XRD = int(data XRD[data XRD.apply(lambda row: row.astype(str).str.contains('No. of points')).any(axis=1)]['Intensity'])
```






## Electrochemical testing

Repository: <a href="mailto:github.com/FedeDat/SuPERCO2">github.com/FedeDat/SuPERCO2</a>

item: data/CuO/electrochemistry

### **Open Circuit Potential**


### **Post-processing**



```
# Read the content of the file
url_OCV = "https://raw.githubusercontent.com/FedeDat/SuPERCO2/main/data/24-10-30_CuO/EC/OCV.csv"
data OCV = pd.read csv(url OCV, header=None, names=['time', 'V'])
pH OCV=float(data OCV['V'][0])
ref_OCV=str(data_OCV['V'][1])
data OCV=data OCV[4:]
data OCV=data OCV.astype(float)
data_OCV['V']=AgCl_to_RHE(data_OCV['V'],6.8,ref_OCV)
plt.figure(figsize=(cm_to_inch(12), cm_to_inch(8)))
plt.plot(data_OCV['time'],data_OCV['V'], )
plt.title('Open circuit potential')
plt.xlabel(r"$t$ (min)")
plt.ylabel(r"$V$ (vs RHE)")
plt.ylim(min(np.min(data_OCV['V']),0.5*np.mean(data_OCV['V'])),max(np.max(data_OCV['V']),1.5*np.mean(data_OCV['V'])))
plt.savefig('OCV.png', format='png', dpi=300, transparent=True, bbox_inches='tight')
#plt.close()
```

## Cyclic voltammetry

### **Post-processing**



### Data analysis & Visualization

```
data_CV_CO2['V']=AgCl_to_RHE(data_CV_CO2['V'],6.8,ref_CV_CO2)
url_CV_N2 = "https://raw.githubusercontent.com/FedeDat/SuPERCO2/main/data/24-10-30_CuO/EC/CV-N2.csv"
data_CV_N2 = pd.read_csv(url_CV_N2, header=None, names=['V', 'j'])
pH_CV_N2=float(data_CV_N2['j'][0])
ref_CV_N2=str(data_CV_N2['j'][1])
data CV N2=data CV N2[4:]
data_CV_N2=data_CV_N2.astype(float)
data_CV_N2['V']=AgCl_to_RHE(data_CV_N2['V'],6.8,ref_CV_N2)
plt.figure(figsize=(cm_to_inch(12), cm_to_inch(12)))
plt.plot(data_CV_CO2['V'],data_CV_CO2['j'], label=r"CO$_2$ saturated",color='red')
plt.plot(data_CV_N2['V'],data_CV_N2['j'], label=r"N$_2$ saturated",color='gray')
plt.title('Cyclic voltammetry')
plt.legend(loc='best')
plt.xlabel(r"$V$ (vs RHE)")
plt.ylabel(r"$j$ (mA/cm$^2$)")
plt.savefig('CV.png', format='png', dpi=300, transparent=True, bbox inches='tight')
#plt.close()
```


#### alculation of CO<sub>2</sub> reduction activity and selectivity

```
# Analyse Faradaic Efficiency
```



## Chronoamperometry

### Post-processing (CA)



### Post-processing (GC)





Work in Products
Work in Products
For liquid Products **Activity & Selectivity** 


**Data analysis & Visualization** 

### Calculation of CO<sub>2</sub> reduction activity and selectivity

```
# Analyse Faradaic Efficiency
  # GitHub API URL for the folder
  api url = "https://api.github.com/repos/FedeDat/SuPERCO2/contents/data/24-10-30 CuO/EC"
  # Get the file listing
  response = requests.get(api url)
  files = response.json()
  metadata_EC = pd.DataFrame(columns=['pH', 'ref', 'V'])
  anl EC = pd.DataFrame(columns=['V RHE', 'Flow Rate', 'Mode', 'Average j', 'Q (mC)'])
  data = pd.DataFrame()
  n1=0
  plt.figure(figsize=(cm_to_inch(18), cm_to_inch(12)))
```

## Chronopotentiometry

### Post-processing (CP)



### Post-processing (GC, HPLC)

1	H2	02	N2	CH4
2	33107.7722			131.4632
3	43991.6529			141.6516
4	52639.5361			148.6601
5	60080.1305			143.4467
6	66500.9825			141.7904
7	71986.6802			139.0178
8	71553.3711			0

1	НСОО-	AcO-
2	3.6034	0.0714
3		
4		
5		
6		
7		



### The FAIR CO<sub>2</sub>R problem

### The SuPERCO2 concept

Open data is limited in CO<sub>2</sub>R

Early data sharing

Zenodo limits interoperability

Characterization: EDX XRD

Al mining on published data

eadm Fle

Electrochemistry: OCV CV CA

CP

images

Add files via up

LICENSE

nitial commit

3 weeks ago

README.md

Next steps

Provide a tutorial

Find sparring partner



## Let's try the power of SuPERCO2!



## Acknowledgements Inspiration Financial support



M.Eng. Moisés Àlvarez (ICIQ)

### Exp collaboration

M.Sc. Ali Zarei (catalyst synthesis, electrochemical testing)

Prof. Simelys Hernández (lab Access)









SuPERCO2 has received funding from the European Union's Horizon Europe research and innovation programme under the Marie Sklodowska-Curie grant agreement No 101104004.

## Thank you! Questions or feedbacks?

## Here or via email (federico.dattila@polito.it)

